



HEMOSURE[®] iFOB TEST

-Baaritaanka Mulikyuulka Kimikada Dhiigga ee Ku Jiro Saxaarada (FIT)



Caafimaad ahaan iyo dhaqaale ahaan waa habka ugu fiican ee lagu Baaro Dhiigga Ku Jiro Saxaarada

"ACG wuxuu taageeraa talada wadajir ahaan la dejiyey oo sheegaayo in habka horey ee baaritaanka saxaarada guaiac la joojiyo iyo loo isticmaalin baaritaanka kansarka (CRC)... ACG waxay ku talinayaan in habka FIT loo isticmaalo baaritaanka kansarka."

Dib loo daabacey kadib markii fasax laga heley MacMillan Publishers Ltd. The American Journal of Gastroenterology (2320970223955), xuquuq daabacad (2009)

Soo koobid & Is-barbardhigid

	HEMOSURE [®] FDA FASAXEY & LAGA DHAAFEY CLIA	BAARITAANKA GUAIIAC EE CAADIGA AH
Habka baaritaanka	Baaritaanka mulikyuulka	Baaritaan Kimiko
Fasiraadda	Natiijo Si Cad loo Akhrin Karo Sida Baaritaanka Urka	Midab baluug, badanaa dhib ah in la fasiro/garto
Cajimaanada Hooseeye ee GI	>96%, Si Sare Gaar u Ah Unugyada Dhiigga Cas ee Insaanka	Hooseeyo, Natiijo qalad ah ayuu ka sheegaa cunto iyo daawo
Dareenka Guud Ahaan	>87%, Lagu ogaado ilaa 50ng hHB/mL ee Unugyada Dhiigga Cas ee Insaanka	50%. Lagu ogaado ilaa 90,000ng hHB/mL ama ka sareeyo oo aan gaar u eheen unugyada dhiigga cas
Saxnaanta	>97%	<86%
Xadayn Daawo Ama Cunto Loo Baahan Yahay Baaritaanka Ka Hor	MIDNA	HAA. Dhowr malmood.
Saambalka loo Baahan Yahay	HAL ayaa loo baahan yahay	SEDDAX ayaa loo baahan yahay



HEMOSURE® iFOB TEST

-Baaritaanka Mulikyuulka Kimikada Dhiigga ee Ku Jiro Saxaarada (FIT)




Ogow: Ma jiraan Xadayn Daawo ama Cunto Baaritaanka Ka Hor

Tilmaamaha Bukaanka

Tallaabada 1 Diyaarinta Saambalka

A Kor u qaad qaybta lagu fariisto musqusha iyo warqadda lagu shubo saambalka ku dheji qarka musqusha. Qaybaha koollada leh si fiican ku dheji geesaha qarka musqusha. Hoos u deji qaybta lagu fashiisto ee musqusha.




B Ku xaar warqadda saambalka dusheeda.



OGOW!

Faah-faahintaada shakhsiyadeed ku qor xaashida ku dhegan tuubada saambalka.




Tallaabada 2 Aruurinta Saambalka


A Fur daboolka baluug-gaduudka ee tuubada saambalka.
DAREERAHA KU JIRO TUUBADA HA KA DAADININ.



B Afka dhuuban dhax-geli ama ku taabo 6 meelo kala duwan ee saxaarada. Hubi inaad caarada sare oo kaliya ku taabato saxaarada. **TUUBADA HA KA BUUXININ, HA KU ARUURIN AMA HA KU FANDHAALIN SAXAARO.**




C Afka dhuuban ee daboolka leh dib u geli tuubada iyo si xoog ah u xir.



Tallaabada 3 Doo Celinta Saambalka

A Soo buuxi boqshadda kuna qor cinwaanka takhtarkaaga ama shaybaarka.



B Tuubada sambalka geli bacda saambalka iyo xir. Bacda saambalka dhax-geli boqshadda iyo xir.



C Rugta caafimaad ama shaybarka si deg-deg ah hab boosto ama laftigaaga keen boqshadda saambalka.



HEMOSURE, INC

Astaanta Hemosure waa mid u diiwaangashan Hemosure, Inc.



DGM36087-1.0